

Diffractometer in Sardana

Teresa Núñez
DESY Photon Science

- Current status
- Details of implementation
- Initialization and use
- Dedicated macros
- Dedicated GUIs
- Next steps



Sardana Workshop
ELI beamlines, 05-06-18

Current status

- Complete diffractometer control from Sardana
- Used at three Petra beamlines (DESY)
- Dedicated GUIs developed in PyQt using Taurus
- Dedicated SEP with status description (not documentation)



Details of the implementation

Diffractometer as a controller device of type PseudoMotor

- Use of the hkl library exclusively done in hkl controller code, not in Sardana core
- Base diffractometer class with all calculations
- Derived diffractometer classes:
 - differing in number of motors and pseudomotors
 - not contains explicit use of hkl library
- Dedicated macros



Details of the implementation (ctd.)

Creating diffractometer controller in Pool:

- Controller type: PseudoMotor
- Controller library: HklPseudoMotorController
- Controller class: depends on number of motors and roles, ex. DiffracE6C
- Motors and Roles: depending on controller class
- Properties:
 - DiffractometerType, depending on geometry, ex. E6C



Details of the implementation (ctd.)

Current diffractometer controller classes:

- Diffrac6C (types: "PETRA3 P09 EH2")
 - motor roles: μ , ω , χ , ϕ , δ , γ
 - pseudomotor roles: h , k , l
- DiffracE6C (types: "E6C", "SOLEIL SIXS MED2+2")
 - motor roles: μ , ω , χ , ϕ , γ , δ
 - pseudomotor roles: h , k , l , ψ , q , α , q_{per} , q_{par}



Details of the implementation (ctd.)

Current diffractometer classes (ctd.):

- DiffracE4C (types: "E4CV", "E4CH", "SOLEIL MARS")
 - motor roles: ω , χ , ϕ , θ
 - pseudomotor roles: h , k , l , ψ , q
- Diffrac4Cp23 (types: "PETRA3 P23 4C")
 - motor roles: ω_t , μ , γ , δ
 - pseudomotor roles: h , k , l , q , α , q_{per} , q_{par} , θ_2 , α_{θ_2} , incidence, emergence



Details of the implementation (ctd.)

Controller device interface:

- Properties:
 - DiffractometerType
- Attributes:
 - parameters describing current diffractometer status
 - commands performing actions, calculations ...
(not extra commands can be defined in controllers)



Sardana / Wiki / SEP4 - Mozilla Firefox

File Edit View History Bookmarks Tools Help

Recibidos - te... hallo - Franzö... 11.13. sqlite3 ... Sardana / ... x Deutschlands... Deutschlands...

sourceforge.net/p/sardana/wiki/SEP4/ Google

- A: a parameter of the current lattice.
- AddCrystal: add a new sample.
- AddReflection: add reflection to current sample.
- AdjustAnglesToReflection: changes the angles associated to the selected reflection.
- AffineCrystal: creates a new sample with '(affine)' attached to the name of the current one and performs the affine. This affine sample is set as the current one.
- Afit: fit value of the a parameter of the current lattice.
- Alpha: alpha parameter of the current lattice.
- AlphaFit: fit value of the alpha parameter of the current lattice.
- B: b parameter of the current lattice.
- Beta: beta parameter of the current lattice.
- BetaFit: fit value of the beta parameter of the current lattice.
- Bfit: fit value of the b parameter of the current lattice.
- C: c parameter of the current lattice.
- Cfit: fit value of the c parameter of the current lattice.
- ComputeTrajectoriesSim: computes the list of trajectories for the current engine and engine mode corresponding to the values of the pseudo axes written in this attribute. The number of arguments has to correspond to the number of pseudo axes for the current engine. The computed trajectories are shown in the TrajectoryList attribute.
- ComputeU: computes UB matrix using the reflections corresponding to the indexes given as arguments.
- Crystal: sample.
- CrystalList: list of samples.
- DeleteCrystal: delete the crystal given in the argument.
- Engine: selected engine. It is taken into account for computing the physical positions corresponding to a movement of a pseudo axis.
- EngineList: list of engines for the diffractometer type corresponding to this controller.
- EngineMode: selected mode for the current.
- EngineModeList: list of the modes corresponding to the current engine.
- HKLModeList: list of the modes corresponding to the hkl engine.
- HKLPseudoMotorList: list of the hkl motor names.
- Gamma: gamma parameter of the current lattice.
- GammaFit: fit value of the gamma parameter of the current lattice.
- LoadReflections: loads reflections for current crystal from ascii file. Reflections will be loaded from the last used file when the Pool is started.
- LatticeReciprocal: reads the values of the reciprocal lattice.
- ModeParametersNames: name of the parameters associated to the current engine mode (if any).
- ModeParametersValues: get/set the value of the parameters associated to the current engine mode (if any).
- MotorList: names of the physical motors associated to the diffractometer.
- ReflectionAngles: angles (computed and theoretical) between the reflections of the current sample.
- ReflectionList: list of reflections for current sample.
- RemoveReflection: remove reflection with given index.
- SaveReflections: saves reflections from current crystal to ascii file. The value written to this attribute is the path to the file, the name of the file is the name of the sample with the termination .ref. If this file already exists a copy will be created adding to the name the current time in seconds.
- SelectedTrajectory: index of the trajectory you want to take when you perform a movement for a given set of pseudo axes positions.
- TrajectoryList: list of trajectories for the current engine and engine mode corresponding to the pseudo axes values written in the ComputeTrajectoriesSim attribute and the engine and engine mode when the calculation was performed. It gives the possibility of checking the trajectories before performing a movement.
- UMatrix: reads current UB matrix values.
- Ux: reads/writes current ux value.
- Uy: reads/writes current uy value.
- Uz: reads/writes current uz value.
- Wavelength.



Initialization and use

- Initialization:
 - sample set to last value or 'default crystal'
 - lattice parameters set to default values
 - geometry created according to DiffractometerType
 - engine set to 'hkl'
- Initial settings can be read from config file:
 - set crystal, ub matrix, reflections, engine and mode
- Movements:
 - using the corresponding motor and role motor devices
- Settings and actions:
 - writing to the controller attributes





Created at 2018-04-29 00:39

DiffractometerType E6C

Crystal srru2o6

Wavelength 4.36871767044

A 5.20795754769 B 5.2079575476 C 5.23299798961

Alpha 90.0 Beta 90.0 Gamma 120.0

R0 0 0.0 0.0 1.0 0 1 0.0 28.90655 91.934 282.357825 -1.52587888991e-08 49.347775

R1 1 0.333333333333 0.0 1.0 0 1 0.0 12.39615 102.452375 282.357825 -1.52587888991e-08 53.17195

Engine hkl

Mode constant_phi_vertical

PsiRef not available in current engine mode

AutoEnergyUpdate 0

U00 0.410 U01 -0.946 U02 0.058

U10 0.067 U11 0.112 U12 1.197

U20 -1.330 U21 -1.016 U22 0.078

Ux -86.2748320555 Uy 2.78981647684 Uz 72.8523304456

SaveDirectory /home/p09user/crystals

Configuration file

Dedicated diffractometer macros

Included in standard sardana macros

- Control diffractometer and display information
- Follow SPEC syntax as close as possible
- Not tested in all diffractometer types
- Examples:
 - br, ubr: move to hkl values
 - hklscan, hscan, kscan, lscan: scans in hkl space
 - setor0, setor1, setorn: set orientation reflections
 - ca, caa, ci: motor positions from hkl values or ψ/χ
 - computeub, setlat, loadcrystal, ...



Dedicated diffractometer GUIs

Implemented using PyQt and Taurus

- Three main GUIs:
 - hkl scans
 - diffractometer alignment
 - UB matrix/lattice parameters
- GUIs connects to:
 - diffractometer controller device
 - MacroServer via Door device





Form



Positions/Limits

H

K

L

Position:

1.44

0.10

0.73

Move to:

1.43

0.10

0.75

Start:

1

1

1

Stop

2

1

1

Angles

omega

chi

phi

theta

2.00

175.92

-147.42

65.29

Parameters

Nb points:

1

Sample time:

1

Select mode:

bisector



Engine:

hkl

Mode:

constant_omega

Start

Stop

Display Angles

Macroserver Connection

Form

Positions/Limits

	H	K	L
Position:	1.44	0.10	0.73
Move to:	1.43	0.10	0.75
Start:	1	1	1
Stop	2	1	1

Angles

omega	chi
2.00	175.92

Parameters

Nb points: 1

Select mode: bissector

Start

Stop

Display Angles

Macroserver Connection

Form <2>

Angles during the scan

omega	chi	phi	theta
3.0000	136.9391	-175.8289	70.5288
3.0000	138.7597	-175.8930	109.4712

Angles

	omega	chi	phi	theta
Position:	2.00	175.92	-147.42	65.29
Move to:	3.00	175.92	-147.42	65.29

hkl

	H	K	L
Position:	1.44	0.10	0.73
Move to:	1.43	0.10	0.75
Select mode:	bissector		
Engine:	hkl		
Mode:	constant_omega		

Scans

Exec Scan:	omega	chi	phi	theta
Range:				
To Max.:	n.n.	n.n.	n.n.	n.n.
Nb points:				

SelectSignal

Stop

Store Refl.

Macroserver Connection

Form

Angles

	omega	chi	phi	theta
Position:	2.00	175.92	-147.42	65.29
Move to:	3.00	175.92	-147.42	65.29

hkl

	H	K	L
Position:	1.44	0.10	0.73
Move to:	1.43	0.10	0.75

Select mode: bissector Engine: hkl Mode: constant_omega

Scans

Exec Scan:	omega	chi	phi	theta
Range:				
To Max.:	n.n.	n.n.		
Nb points:				

Stop Store Refl.

Form <2>

Select:

exp_t01

Signal: exp_t01

Sample Time:

SelectSignal

Macroserver Connection

Herzlich willkommen!

Form

UB Matrix

UB11	4.07999045921	UB12	-2.49827362821e-16	UB13	-2.49827362821e-16
UB21	0.0	UB22	4.07999045921	UB23	-2.49827362821e-16
UB31	0.0	UB32	0.0	UB33	4.07999045921

U Vector

Ux	0.00	Uy	0.00	Uz	0.00
	0.00		0.00		0.00

Lattice

a	1.54	b	1.54	c	1.54
	1.54		1.54		1.54
alpha	90.00	beta	90.00	gamma	90.00
	90.00		90.00		90.00

Update

ComputeU

Reflections List

Edit Reflections

A decorative illustration of a couple dancing under a street lamp at night. The man is wearing a light-colored shirt and dark pants, and the woman is wearing a dark dress. They are both smiling and looking at each other. The street lamp is tall and has a warm glow. The background is dark with some light streaks.

Herzlich willkommen!

Form

UB Matrix

UB11	4.07999045921	UB12	-2.49827362821e-16	UB13	-2.49827362821e-16
UB21	0.0	UB22	4.07999045921	UB23	-2.49827362821e-16
UB31	0.0	UB32	0.0	UB33	4.07999045921

Form <2>

	Index	H	K	L	omega	chi	phi	tth
Ux	0	1.0000	1.0000	1.0000	30.0000	0.0000	87.1000	22.0000
	1	1.0000	0.0000	1.0000	28.0000	1.0000	86.5000	23.1000
	2	1.0000	0.0000	0.0000	45.0000	23.0000	4.5000	10.0000
Lattice								
a								
alpha								

Update

Clear Apply

Next steps

- Test of current implementation:
 - performance
 - completeness
- Discuss about other possible implementations

