



# Continuous Scans in Sardana (SEP6)

Teresa Nunez,  
Zbigniew Reszela -  
on behalf of the Alba Controls Group  
Sardana Workshop 2015



*Motion & acquisition during the step scan.*



*Motion & acquisition during the continuous scan*

# Important assumptions

- All elements, but slave motors, must be defined in the same Pool.
- Support only linear trajectories – constant velocity.
- Global integration time per scan point record.

# Important assumptions

- All elements, but slave motors, must be defined in the same Pool.
- Support only linear trajectories – constant velocity.
- Global integration time per scan point record.
- Software development approach:
  - reuse the most from the existing Sardana
  - iterative development instead of upfront plan and designs
  - let's first obtain the complete model and later optimize it
  - strong emphasis on automated tests

Only one motor is used as position synchronization source.

Continuous velocity move until the last sample is acquired

The same measurement group must be reusable between step and continuous scans (if hardware allows that).


`ascanc <motor> <start_pos> <final_pos> <nr_interv> <integ_time>`


#Pt	No	motor	integ_time	expchannel	dt
0		13.0	0.01	8.780	0.0
1		13.1	0.01	8.698	0.2


Nominal motor positions


NEW:  
integration time  
set point

Nominal delta time  
from the start












Hardware	Supported domains	Non equidistant configuration	Config. active interval – gate support	Interpret configuration as absolute	Start on position and cont. in time	Require position calibration
<b>Software</b>	Time Position	YES	YES Time Position	YES	YES	<b>NO</b>
<b>Icepap</b>	Position	YES	<b>NO</b>	YES	NO	YES
<b>Zebra</b>	Time Position	NO	YES Time Position	YES (Position)	YES	YES (Internal)
<b>TurboPmac2 (standard feature)</b>	Position	NO	YES Position	YES	NO	YES
<b>TurboPmac2 (PLC)</b>	Time Position	YES	YES Time Position	YES	YES	?
<b>NI6602</b>	Time Position*	NO	YES	<b>NO</b>	NO*	YES
<b>TFG</b>	Time	YES	YES	<b>NO</b>	-	-




*Group*

```
[{Delay: {Time: 0.3, Position: 400},
  Initial: {Time: None, Position: 0},
  Active: {Time: 0.1, Position: 10},
  Total: {Time: 0.15, Position: 15},
  Repeats: 1000},
...]
```

- Continuous scan fills the synchronization parameters in position and time domains.
- Time scan and step scan\* fills just the time domain.
- Step scan uses just one *group* and starts the measurement group multiple times\*.
- TriggerGate controller chooses the most appropriate parameters – controlled by users with extra attributes.
- Synchronization with first trigger on position and continues in time are possible.
- General rule: position domain takes precedence for Initial and Total but time domain takes precedence for Active.

\* step scans could use synchronization with Repeats > 1 in the future



- Useful in software synchronized channels – helps to avoid skipped acquisitions
- $\text{wait time} = \max(\text{re-arm\#1}, \text{re-arm\#2}, \text{MG wait time})$
- Affects: motors velocities, total interval (time)

# Acquisition – hardware

Hardware	Type of controller	Multiple channel	Channels are independent	Synchroni- zation	Require external timer	Allows readouts while acquiring
<b>Software (TaurusAttr.)</b>	Counter, 0D, 1D, 2D	YES	<b>YES</b>	SW Trigger SW Gate	NO	YES
<b>AlbaEM</b>	0D (Electrometer)	YES	NO	HW, SW Trig. SW. Gate	NO	<b>NO</b>
<b>Adlink2005</b>	0D (ADC)	YES	NO	HW, SW Trig. HW, SW Gate	NO	YES
<b>Keithley</b>	0D (Electrometer)	YES	NO	HW, SW Trig. ?, SW. Gate	NO	<b>NO</b>
<b>NI6602</b>	Counter, Timer	YES	<b>YES</b>	HW, SW Trig. HW, SW Gate	<b>YES (can use internal)</b>	YES
<b>Mythen</b>	1D	YES (but used as single)	NO	HW, SW Trig. HW, SW Gate	NO	YES
<b>CCDs &amp; Detectors</b>	2D	<b>NO</b>	NO	HW, SW Trig. HW, SW Gate	NO	YES
<b>SIS3820</b>	1D	YES	NO	HW, SW Trig.	<b>Yes</b>	<b>NO</b>
<b>XIA</b>	1D	YES	NO	HW, SW Trig.	NO	<b>NO</b>
<b>Zebra</b>	Counter	YES	NO	HW Trig.	NO	YES (but not in use)


- Start action:
  - configure synchronization: software trigger, software gate, hardware trigger or hardware gate
  - configure number of repetitions
  - call the Load sequence on all triggered channels
  - call the Start sequence in the reversed order of configuration
- Action loop:
  - while channels are acquiring
    - call the State sequence
    - call the Read sequence (for channels with the online readout capability)
  - call the Read sequence (for channels without the online readout capability)

- SardanaValue contains: value, timestamp and **index**.
- Controllers may return SardanaValue objects or just the python base objects: float, sequence of floats, or numpy objects (in the last case the acquisition action fills the rest of the fields).
- GSF receives data in chunks and fills the records in base of the indexes.
- Zero order interpolation is applied in case of missing data.
- Interpolated data must be easily distinguishable from the raw data.

# Configuration up to now

Some controllers create Timer just to pass the integration time to the HW

Timer, Monitor, TGElement are competing for the synchronizer role



Channel	enabled	output	Shape	Data Type	Plot Type	Plot Axes	Timer	Monitor	Trigger	Conditioning	Normalization	NeXus Path	Synchronization
nictpost	true		[]	float64	Spectrum	<mov>	nictpost	nictpos01	Trigger		No		itg01
nictpos01	true	true	[]	float64	Spectrum	<mov>	nictpost	nictpos01	Trigger		No		itg01


- Timer & Monitor, master Timer & Monitor roles are not clear to all users
- Many times Timer/Monitor element do not change between the measurement groups
- expconf does not allow to use Timer or Monitor from another controller
- expconf forces the same Timer and Monitor for the whole controller
- expconf forces the same synchronization (trigger or gate) for the whole controller



# Config. w/o Timers/Monitors

Some controllers create Timer just to pass the integration time to the HW

Timer, Monitor, TGElement are competing for the synchronizer role



Channel	enabled	output	Shape	Data Type	Plot Type	Plot Axes	Timer	Monitor	Trigger	Conditioning	Normalization	NeXus Path	Synchronization
nictpost	true	true	[]	float64	Spectrum	<mov>	nictpost	nictpos01	Trigger		No		itg01
nictpos01	true	true	[]	float64	Spectrum	<mov>	nictpost	nictpos01	Trigger		No		itg01


Less Timers means:

- less elements to control in the acquisition actions
- less Tango devices & events
- no problem with the Value attribute format
- less columns in the records

Configuration w/o Timer/Monitor means:

- simpler expconf GUI
- simpler data structures for the programmer

- Configuration is a table containing **all** information about the acquisition and synchronization elements
- **Timer** and **Monitor** columns should disappear in favor of the **Synchronizer** column
- No Timer for the ExpChannels means no info about the real integration time, is it a problem?
- The Synchronizer column represents the element in charge of controlling the acquisition (either software or hardware)
- Valid options for the Synchronizer column are any of: **TriggerGate** elements or **CounterTimer** elements



Channel	Control	Synchronizer
CounterTimer#1	Trigger	HW TriggerGate
ExpChannel#1	Trigger	HW TriggerGate
CounterTimer#2	Gate	CounterTimer#1
ExpChannel#2	Gate	CounterTimer#1
ExpChannel#3	Trigger	SW TriggerGate

Timers **do not** belong to the synchronization action

Timer


Slow channel

Call **LoadOne**  
on all of them

Direction of the  
synchronization  
control



**Icepap**



Channel	Control	Synchronizer
ExpChannel	Gate	CounterTimer
CounterTimer	Trigger	HW Position TG

