

Sardana at Petra III (DESY)

Teresa Núñez
DESY Photon Science

- Current status
- Installation
- Sardana implementation
- Configuration and Startup
- Implementation specific procedures
- Taurus GUIs
- Extra features



Sardana Workshop
DESY, 06-09-15

Current Status

- Sardana installed in all DESY-PETRA beamlines
- Used:
 - primordially at Dynamics, Extreme Conditions, Resonant Scattering and Diffraction, Coherence Applications
 - partially at Variable Polarization XUV, Hard X-ray Micro/Nano-Probe, High Resolution Power Diffraction, HighRes Diffraction
 - none at MINAXS, High Energy Material Science



Installation

Own debian packages merging Sardana develop branch with Sardana-DESY repository with own modifications

Sardana develop



Sardana-DESY:


- nxs/fio recorders
- General hooks
- General on_stop
- Check-condition
- Handling motor limits



Debian package



Sardana Implementation




Restricted by naming convention and MG


Configuration and Startup

- Configuration based on xml file:
 - all devices
 - measurement groups
 - diffractometer with initial configuration
- Dedicated tools for starting/stopping Sardana and performing any user required action at startup.

Sardana is started/configured from scratch (related DB info removed and created again) if changes in the xml file



Configuration and Startup



Implementation specific procedures

Macros running from Spock

- DESY specific macros:
 - energy macros: `escan`, `me`, `escanexafs`, ...
 - hkl macros: `hklscan`, `br`, `wh`, `setorn`, `luppsi`, ...
 - petra macros: `wait_for_petra`, `wait_for_beam`, ...
 - motor limits macros: `hasy_adjust_lim`, ...
- Beamline specific macros:
 - Maia macros
 - Zebra & XIA macros
- User specific macros



Taurus GUIs

- General GUIs occasionally used:
 - Macrogui (expconf)
 - imagegui
- Beamline specific GUIs:
 - motor display
 - full beamline control (Variable polarization XUV)



Exp2 beamline panel

overview

RMU2 control

status

status 2

screens

vacuum interlock

beamshutter

PETRA III info

P04 elog

PETRA III elog

Help

PETRA III Current

100.4843mA

Optic ready

screens out

ps out

scan active

info field

('No Scan active...')

Emergency Stop

Master energy

(Mono and Undulator)

+

1

0

0

0

0

0

0

0

0

go

Mono energy

779.295eV

+

0

7

7

9

3

0

0

mono Cff

2.079

Mirror Name

PlaMir Pt

Grating Name

G1200nm12

UNDUALTOR Energy

777.98eV

+

0

7

7

8

0

0

Desired Gap

22225μm

+

2

2

2

2

5

0

0

go

UNDUALTOR GAP

22225μm

UNDUALTOR Shift

-17680μrad

The undulator has reached the requested position and is waiting for a new request.

Beamline Motor

Exit Slit

48

+

0

0

0

0

4

8

micron

exsu2bpm

-1.00

-

0

0

1

0

0

mm

exsu2baffle

30.01

+

0

3

0

0

1

mm

Reset

Apply

EXSU2 motor breaks

trans break

off

on

slit break

off

on

jlive

print motor pos to elog

show astor exp1

load config

show astor exp2

save config with name

show panel terminal

Keithley panel

archive viewer

scan panel

password expert panel


open expert panel



Beamline Control with Taurus (p04- F. Scholz)

Herzlich willkommen!






Beamline Control with Taurus (p04- F. Scholz)

Herzlich willkommen!






Beamline Control with Taurus (p04- F. Scholz)

Herzlich willkommen!






Beamline Control with Taurus (p04- F. Scholz)

Herzlich willkommen!





Beamline Control with Taurus (p04- F. Scholz)

Extra features

Requested and in use at DESY

- General scan hooks
- General macro abort function
- General condition for repeating scan points

Add flexibility to scans without defining new macros

